



May 1, 2009

"Going beyond expectations"


DEALER INFORMATION

BILL TO: (For Invoicing, if Different then "Ship To" Location)		SHIP TO:	
Dealer Name:	_____	Dealer Name:	_____
Contact:	_____	Contact:	_____
Address:	_____	Address:	_____
City:	_____	City:	_____
Province:	_____	Postal Code:	_____
Phone Number:	_____	Phone Number:	_____
Fax Number:	_____	Fax Number:	_____
Email Address:	_____	Email Address:	_____
Terms:	_____	Secondary Discount:	0%
Primary Discount:	0%	Total Discount:	0.00%

ORDER / QUOTE INFORMATION

Date: _____ Quote/Order Number: _____
P.O. Number: _____ Account Number: _____

INFORMATION ON THE USER

LEFT	RIGHT		
A _____	_____		User name or number: _____
B _____	_____		User weight: _____
C _____	_____		User height: _____
D _____	_____		Seat cushion thickness: _____
E _____	_____		
F _____	_____		
G _____	_____		
H _____	_____		
J _____	_____		
K _____	_____		

POWER BASE SELECTION

Wheelchair Model	¹ STFH - Tilt & Lift	¹ STFH - Tilt
<input type="checkbox"/> Quickie Xperience	19" (B & E)	16 3/4" (A)
<input type="checkbox"/> Quickie Xplore	19 1/4" (B & E)	17" (A)
<input type="checkbox"/> Quickie Rhythm	19 1/4" (A & E)	17 1/4" (A)
<input type="checkbox"/> Quickie Groove	19 1/2" (A & E)	17 1/2" (A)
<input type="checkbox"/> ² Quickie Freestyle F11	20 1/4" (A & E)	18 1/4" (B)
<input type="checkbox"/> ² Quickie Freestyle M11	19 1/4" (A & E)	17 1/4" (B)
<input type="checkbox"/> Quickie S646 SE <input type="checkbox"/> Std. <input type="checkbox"/> Wide	20 1/2" (B & E)	18 1/4" (C)
<input type="checkbox"/> ³ Quickie S525 <input type="checkbox"/> Narrow <input type="checkbox"/> Std. <input type="checkbox"/> Wide	N/A	18" (C)
<input type="checkbox"/> ³ Quickie P220 / P222 <input type="checkbox"/> Std. <input type="checkbox"/> Wide	N/A	19 1/2" (D)
<input type="checkbox"/> Invacare TDX SP (Group 22 Batteries)	20 1/2" (A & E)	18" (B)
<input type="checkbox"/> Invacare TDX SP (Group 24 Batteries)	20" (A & E)	17 1/2" (B)

Wheelchair Model	¹ STFH - Tilt & Lift	¹ STFH - Tilt
<input type="checkbox"/> Invacare Storm 3G <input type="checkbox"/> Short <input type="checkbox"/> Long	20 1/2" (B & E)	18" (C)
<input type="checkbox"/> Invacare TDX SR	20 1/2" (A & E)	18" (B)
<input type="checkbox"/> Invacare TDX SC (Spree)	TBA (A & E)	TBA (B)
<input type="checkbox"/> Invacare <input type="checkbox"/> TDX3 <input type="checkbox"/> TDX4 <input type="checkbox"/> TDX5	20" (A & E)	17 3/4" (A)
<input type="checkbox"/> Invacare Pronto M91	20" (A & E)	18" (B)
<input type="checkbox"/> Quantum 610	20" (A & E)	18" (B)
<input type="checkbox"/> Quantum <input type="checkbox"/> 600E <input type="checkbox"/> 600 <input type="checkbox"/> 600XL	20 1/4" (A & E)	18" (A)
<input type="checkbox"/> Quantum <input type="checkbox"/> 6000 <input type="checkbox"/> 6000Z <input type="checkbox"/> 6000XL	19 3/4" (A & E)	17 1/4" (B)
<input type="checkbox"/> Quantum R4000	20" (B & E)	17 1/2" (A)
<input type="checkbox"/> Quantum 1121	20 1/2" (A & E)	18 1/2" (B)
<input type="checkbox"/> Quantum 1122	20 1/2" (A & E)	18 1/2" (B)

FOR INTERNAL USE ONLY
Note: For a tilt only with a Large seat depth range between 18" and 22", please replace the C-5200 tilt and use the C-6400 tilt. Seat to floor height may increase by up to 3/4".
A - C6400 C - C6000 E - C4200
B - C5200 D - C6100

Other Power Wheelchair Base. Specify: _____

ELECTRONICS (INCLUDED WITH POWER WHEELCHAIR)

QUICKIE	PRIDE - QUANTUM	INVACARE - MODULE	JOYSTICK
<input type="checkbox"/> Qtronix	⁴ <input type="checkbox"/> Q-Logic	<input type="checkbox"/> MKIV-A	
<input type="checkbox"/> Delphi QR	<input type="checkbox"/> Q-Logic NE	<input type="checkbox"/> MKIV-RII	
<input type="checkbox"/> Delphi QC	<input type="checkbox"/> Remote Plus	<input type="checkbox"/> MK5-EX	<input type="checkbox"/> MPJ
<input type="checkbox"/> Standard VSI	<input type="checkbox"/> VR2	<input type="checkbox"/> MK5-NX	<input type="checkbox"/> SPJ
<input type="checkbox"/> 6-Key VSI	<input type="checkbox"/> 6-Key VR2	<input type="checkbox"/> MK690ACC	
<input type="checkbox"/> R-net	<input type="checkbox"/> Standard VSI	<input type="checkbox"/> MK6TT <input type="checkbox"/> GTRAC	
<input type="checkbox"/> 6-Key VR2	<input type="checkbox"/> 6-Key VSI	⁵ <input type="checkbox"/> MK690	

FOOTNOTES

- (STFH) Seat to floor height was measured at the front of the seat pan with a 0° Pre-Tilt with the following seats: Invacare: ADJASBA seat ----- Pride: Synergy seat ----- Quickie: Std. rehab
- For the Quickie Freestyle F11 & M11 with the standard rehab seat, cutting part of seat is required. See install manual for directions and details.
- For the Quickie S-525 and the P220 / P222 with VSI electronics, they must be tilt ready and equipped with a female Redel "Pigtail" connector in order to obtain 24V power and drive lock-out.
- When a chair is equipped with an 11" seat lift and Q-Logic electronics, reduced speed requires going through drive controls. Please select one of these part numbers for reduced speed: SWTJSS or SWTJSD.
- MK690 is not compatible to be used with power seating options.



INSTALLATION		Retail Price	Discount Price	Net Price
<input type="checkbox"/>	NOTE-IN-HOUSE	Factory Installation (The Wheelchair will drop-ship from the manufacturer)	STD	
<input type="checkbox"/>	NOTE-KIT	Dealer - Ready to Install Kit	N/C	

INTERFACE POWER MODULES				
<input type="checkbox"/>	PTS50CG-	6	50° CG Power Tilt Interface	\$3,495
<input type="checkbox"/>	PES-C4200	6, 7	11" ADD-ON Power Seat Lift Interface	\$3,495

ELECTRONICS - SINGLE (1) OR DUAL (2) ACTUATOR CONTROLLERS				
<input type="checkbox"/>	SAC	Single Actuator Controller (Includes Reduced Speed with the Power Elevating Seat)		STD
<input type="checkbox"/>	DAC	Dual Actuator Controller (Includes Reduced Speed with the Power Elevating Seat)		\$495

Please select one activation and position from the list below

<input type="checkbox"/>	SWTOGGLE	8 or 9	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC) - <input type="checkbox"/> User <input type="checkbox"/> Attendant	STD
<input type="checkbox"/>	SWPUSH1	8	Single Push Button Switch (Single Mode: SAC) - <input type="checkbox"/> User <input type="checkbox"/> Attendant	N/C
<input type="checkbox"/>	SWPUSH2	8 or 9	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC) - <input type="checkbox"/> User <input type="checkbox"/> Attendant	N/C
<input type="checkbox"/>	SWTJSS-	8	Single Actuator Activation through Wheelchair Controls	\$995
<input type="checkbox"/>	SWTJSD-	9	Dual Actuator Activation through Wheelchair Controls	\$995
<input type="checkbox"/>	SWTVR2-	10	Tilt only activation through 6-key P&G VR2 joysticks	N/C
<input type="checkbox"/>	SWTVSI-	11	Tilt only activation through 6-key P&G VSI joysticks	N/C
<input type="checkbox"/>	SWTMK6-ACCT	12	Tilt only activation through Invacare: <input type="checkbox"/> SPJ-ACC <input type="checkbox"/> SPJ-PSS ¹³ <input type="checkbox"/> MPJ ¹⁴	N/C

JOYSTICK MOUNTING LOCATION				
<input type="checkbox"/>	Left hand side		<input type="checkbox"/>	Right hand side

SEAT SELECTION						
<input type="checkbox"/>	Quickie Standard rehab seat	<input type="checkbox"/>	Quickie A.S.A.P. seat	N/C		
<input type="checkbox"/>	Quickie Perfect Fit	<input type="checkbox"/>	Quickie Eclipse Seat	N/C		
<input type="checkbox"/>	Pride Synergy seat	<input type="checkbox"/>	Pride Specialty seat	N/C		
<input type="checkbox"/>	Invacare ADJASBA - Adjustable Pediatric (12" to 16" wide) (12" to 18" deep)			N/C		
<input type="checkbox"/>	Invacare ADJASBA - Adjustable Adult (16" to 20" wide) (16" to 22" deep)			N/C		
<input type="checkbox"/>	Invacare ADJASBA - Adjustable Adult (20" to 24" wide) (16" to 22" deep)			N/C		
<input type="checkbox"/>	Invacare ASBA - Adjustable Pediatric (12" to 16" width & depth) / Adult (14" & 15" wide) (16" to 22" deep)			N/C		
<input type="checkbox"/>	Invacare ASBA - Non - Adjustable Adult (16" to 24" wide) (16" to 22" deep)			N/C		
<input type="checkbox"/>	¹⁵ Pride Captain seat	<input type="checkbox"/>	¹⁵ Invacare Captain seat	<input type="checkbox"/>	¹⁵ Quickie Captain seat	N/C

Specify seat width: " Specify seat depth: "

FRONT RIGGINGS				
<input type="checkbox"/>	SOMFR	16	Using the OEM footrests from: <input type="checkbox"/> Invacare <input type="checkbox"/> Pride <input type="checkbox"/> Quickie	STD
<input type="checkbox"/>	SOAAF90AO	17	90° Fixed Center Mount Interface Bracket for the Flip-up one piece Foot Platform	\$225
<input type="checkbox"/>	SOFUFP		Flip-up one piece Foot Platform	\$175

TOTALS				Retail Price	Discount Price	Net Price
				\$0.00	\$0.00	\$0.00

NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS

FOOTNOTES:

- 6 The maximum user weight capacity is 300 Lbs.
The 11" power seat lift comes standard with reduced speed set at 30% of the preset speeds.
- 7 **Notice on Pride Q-Logic electronics: Reduced speed requires going through drive controls**
Please select one of these part numbers for reduced speed: SWTJSS or SWTJSD.
- 8 Requires the Single Actuator Controller (SAC)
- 9 Requires the Dual Actuator Controller (DAC)
- 10 Requires the P&G 6-Key VR2 Joystick from Quickie or Pride.
- 11 Requires the P&G 6-Key VSI Joystick from Quickie or Pride.
- 12 Requires the SPJ-ACC, SPJ-PSS or MPJ Joystick with the MK690-ACC controller from Invacare.
- 13 Requires the use of a Single Switch with a 1/8" phono plug (Mono)
- 14 Requires the use of a Single or a Dual Switch with a 1/8" phono plug
- 15 The 50° CG power tilt interface option is not available with a captain seat. In-House installation required with all captain seats
- 16 Amysystems does not supply OEM footrests from either Quickie, Invacare or Pride.
- 17 Allows the footboard to follow the seat elevator and/or tilt. This is not a standard option for the following power chairs:
Storm 3G, S646, S525 or R-4000. Depending on the configuration, the option may be available as a custom.

Amysystems
1650 Chicoine
Vaudreuil-Dorion, QC.
J7V 8P2

Phone: (888) 453-0311
Fax: (877) 501-8458

service@amysystems.com
www.amysystems.com
Page 2 of 2